



**soniceraser<sup>™</sup> pro**  
by MICHAEL TODD BEAUTY

World's First  
**SONIC INFUSER**  
With Triplex Infusion Technology<sup>™</sup>  
Sonic | Ionic | Thermal  
Face, neck and décolletage

**INSTRUCTIONS  
FOR USE**

MICHAEL TODD



# CONGRATULATIONS

Dear Customer,

Thank you and congratulations on purchasing **soniceraser pro**, a revolutionary infusing tool that makes your favorite skin care brands work even better and faster. Using our unique Triplex Infusion Technology,<sup>™</sup> the **soniceraser pro** combines 3 technologies, working together, to improve the absorption of the active ingredients in your skin care products up to **8x better** than with hands alone.

With **soniceraser pro**, thermal technology gently warms the skin to open pores and help prepare the skin to better receive and absorb the active ingredients in your skin care products.

Next, sonic technology operating at 36000 *movements per minute* quickly, painlessly and safely helps push your skin care products into the skin, resulting in increased absorption.

Finally, Ionic technology in the **soniceraser pro** eliminates the need for an On/Off switch. The Ionic technology forms a circuit when touched to your skin to conveniently turn the **soniceraser pro** "on" and then "off" when removed from the skin avoiding splatter or wasting product.

Enjoy younger looking skin of the face and décolletage with anti-aging **soniceraser pro** with unique Triplex Infusion Technology<sup>™</sup>.

*For best results, read these instructions for use completely before using the **soniceraser pro** for the first time. If you have any questions or problems with your device, contact us at Michael Todd Beauty Customer Service at : 1.772.800.7654 or [customerservice@michaeltoddbeauty.com](mailto:customerservice@michaeltoddbeauty.com)*

## REGULATORY INFORMATION

U.S. This Class B digital product complies with Part 15 of the U.S. Federal Communications Commission (FCC) rules. Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Canada. This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference and (2) this device must accept any interference including interference that may cause undesired operation of the device.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes: (1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement."

Mfg: Michael Todd Beauty

Model: **soniceraser** *pro*

Device complies with  
CAN ICES-3 (B)/NMB-3(B)

For home or office use. Not intended for use in machinery, medical or industrial applications.

### Radio and TV Interference Regulations.

This product has been tested and found to comply with the limits for a Class B digital device pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This product generates, uses and can radiate radio frequency (RF) energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications devices (such as televisions, baby monitors, cordless phones, etc). However, there is no guarantee that RF interference will not occur in a particular installation. If this product does cause harmful interference to radio or television reception, which can be determined by turning the product off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antennae of the other radio-communications device
- Move the equipment farther away from the radio or tv
- Connect the product into an outlet on a circuit different from that to which the other device is connected
- Consult the manufacturer's customer service or an experienced radio/tv technician for more suggestions.

# IMPORTANT SAFETY INFORMATION

**BASIC SAFETY PRECAUTIONS.** When using electrical products, especially when children are present, basic safety precautions should always be taken. IF YOU THINK THE POWER IS OFF WHEN THE PRODUCT IS OFF, YOU ARE WRONG. If the plug is in, the power is on. Unplug electrical appliances when you are not using them.

**DANGERS. To reduce the risk of electrocution:**

1. **DO NOT** get the Handle wet or immerse in water. **DO NOT** place or store this product where it can fall or be pulled into a tub or sink.
2. **DO NOT** place Handle, Cradle or Charging Cord in or drop it into water or other liquid. The Handle, Cradle or charging Cord are not for immersion or for use in the shower or bath.
3. **DO NOT** reach for the Handle, Cradle or Charging Cord that has fallen into water. Unplug immediately.
4. **DO NOT** use the Charging Cord where the cord while plugged in can come in contact with water.
5. Keep this product away from open flame or flammable aerosol products. The device is not suitable for use in the presence of a flammable anesthetic mixture with air or oxygen or nitrous oxide.
6. **DO NOT** use this product on broken or abraded skin or if you have skin reaction hypersensitivity.

**WARNINGS. To reduce the risk of burns, electrocution, fire or injury to persons:**

1. Close supervision is necessary when this product is used by, on, or near children or physically challenged persons.
2. Use this product only for its intended use as described in these Instructions for Use. **DO NOT** use attachments other than those recommended and provided by Michael Todd Beauty.
3. Keep the Charging Cord away from heated surfaces.
4. Never use this product while sleeping or drowsy.
5. **DO NOT** use this product where oxygen is being administered.
6. Never operate this product if it is damaged.
7. **DO NOT** use the Charging Cord for recharging products other than this product or if the cord or either end of the cord is damaged, is not working properly, or has been dropped into water.
8. The product has no serviceable parts. Never attempt to open the Handle for any reason. If the Charging Cord is damaged, it must be unplugged and replaced. Contact customer service at 772.800.7654 or email us at [customerservice@michaeltoddbeauty.com](mailto:customerservice@michaeltoddbeauty.com) for instructions on how to obtain a replacement.

9. Use this product with the supplied Charging Cord only.
10. **DO NOT** disassemble, heat above 60° C (140°F) or incinerate.

#### IMPORTANT CAUTIONS.

1. This is not a medical device and is not intended to treat medical conditions. If you are under a doctor's care for any skin condition, if you think you may have a skin condition or if you have very sensitive skin, consult your physician before using this product.
2. This product uses magnets which may affect internal medical devices such as pacemakers. Keep away from internal medical devices. To avoid damage, keep magnetic media such as computer discs, credit cards and tapes away.
3. If you experience any pain or discomfort from using this product, discontinue use and consult a physician.
4. This device is electrical and contains small parts. Keep out of reach of children. THIS PRODUCT IS NOT DESIGNED OR INTENDED PRIMARILY FOR CHILDREN 12 YEARS OF AGE OR YOUNGER.

# CONTENTS

Your soniceraser *pro* comes with the following:

- A. soniceraser *pro* Infusion Tool
- B. US Wall Plug
- C. Charging-Storage Cradle



# GETTING TO KNOW YOUR SONICERASER *pro*



## SONICERASER *pro* FEATURES

Products for the delicate skin of the face and décolletage can work better when absorbed deeper and more evenly into the skin. With our unique Triplex Infusion Technology,<sup>™</sup> **soniceraser *pro*** uses a combination of sonic, ionic and thermal technology to gently, safely and painlessly infuse the ingredients of your skincare products deeper into the skin where they can do the most good.

First: As you use the **soniceraser *pro***, the Thermal technology in the Treatment Head heats to a comfortable 42° C to open pores and help prepare the skin for enhanced absorption of the active ingredients in your skin care products.

Next: the Sonic technology in the **soniceraser *pro*** applies the beneficial ingredients in your anti-aging serums and creams treatment products with gentle micro-massaging pulses up to 36000 movements per minute. Aestheticians recommend that products used on the delicate skin of the face and décolletage be applied with gentle tapping rather than rubbing to avoid irritating or injuring the skin. The gentle sonic micro-massaging pulses of the **soniceraser *pro*** simulate the gentle tapping recommended by aestheticians without irritating rubbing which can injure the skin.

Finally: Ionic technology in the **soniceraser *pro*** eliminates the need for an On/Off switch. The Ionic technology forms a circuit when touched to your skin to conveniently turn the **soniceraser *pro*** "on" and then "off" when removed from the skin avoiding splatter or wasting product.

The **soniceraser pro** is cordless and rechargeable and comes with a large Treatment Head ideal for applying anti-aging serums and creams to the delicate skin of the face, including the areas around the eyes and lips, as well as the décolletage which includes the upper chest area, neck, cleavage and shoulder areas. The **soniceraser pro** will:

- Leave the skin feeling energized and refreshed
- Provide immediate and long-lasting hydration
- Brighten, firm and tighten, and
- Reduce the appearance of fine lines, dark circles and under eye puffiness.

# OPERATING INSTRUCTIONS

1. The soniceraser *pro* is intended for use on the skin of the face and décolletage. It is gentle enough for daily use, morning and night. You should feel a pleasant warm and tingling sensation when using the device. This is normal. Warning: Stop using if discomfort or irritation occurs. To prevent eye injury, do not use or place the Treatment Head directly onto the eye or the upper eyelid.
2. Prior to using the soniceraser *pro* for the first time, fully charge the unit following the instructions contained in this User Guide.
3. The Ionic technology in the soniceraser *pro* turns the device on and off automatically so no switch is required. When you hold the soniceraser *pro* with your hand touching the Sensor the device will turn on when you touch the Treatment Head to your skin and shut off when it is removed.
4. Wash your face following your normal routine. Pat dry. You are ready to begin your treatment. Apply a pearl-sized amount of anti-aging cream or serum directly to the area to be treated or to the Treatment Head.

**Important: If your device continues to run while not in use, please note to clean and wipe any excess product that may be left over on your soniceraser *pro* Treatment Head.**

5. Holding the Handle comfortably in your hand, gently position the Treatment Head on clean, dry skin. Using gentle pressure, glide the Treatment Head in small overlapping circular motions across the area being treated until the serum or cream is absorbed. Complete one area then go on to the next treatment area until finished.
6. Following treatment, clean the Treatment Head with a clean towel or tissue and store in the Storage Cradle provided. Charge as needed.



# CHARGING

1. Universal Charging Cradle. Your *soniceraser pro* comes with a universal Charging Cradle to both hold and recharge the rechargeable Handle. The universal Charging Cradle is compatible with international and domestic outlets (100-240 Volts AC; 50/60 Hz) with the respective country plug adapter. If the shape of the plug does not fit the power outlet, use an attachment plug adapter of the proper configuration for the power outlet.
2. To charge the *soniceraser pro*, plug the Charging Cord into a wall outlet and plug the DC jack into the back of the Charging Cradle. Place the Handle into the Charging Cradle as shown. The Charge Indicator light will light up red indicating the unit is charging and will turn green when the unit is fully charged. The Charging Cord may become warm in normal operation while charging. Unplug the Charging Cord from the wall and the Charging Cradle when not in use. DO NOT unplug the Charging Cord from the wall outlet when it is still connected to the unit as this will cause discharging of the battery.
3. Prior to initially using your *soniceraser pro*, it must be fully charged. Allow 8-10 hours for initial and subsequent recharging. A full charge allows up to 40 minutes of continuous use. DO NOT charge your *soniceraser pro* for more than 40 continuous hours.



4. Infrequent Use. The Handle needs to receive a charge occasionally. If you rarely use the Handle, recharge it in the Charging Cradle at least every 30 days to maintain its charging capability. NOTE: Failure to regularly recharge the battery may cause the battery to stop recharging. We recommend that you keep the Handle in the Charging Cradle when not in use.
5. If your **soniceraser pro** stops working, charge it in the Charging Cradle for a minimum of 5 hours. If the issue is not resolved, please contact Customer Service.

## CARE AND MAINTENANCE

1. **Cleaning.** The Handle and Treatment Head can be cleaned using a clean towel or tissue to remove residue buildup. **DO NOT** get the Handle or Charging Cradle wet or immerse in water. Unplug the Charging Cradle before cleaning it. After unplugging the Charging Cradle, wipe with a soft, damp cloth.
2. **Rechargeable Battery.** Prior to initially using your *soniceraser pro*, it must be fully charged. Allow 8-10 hours for initial and subsequent recharging. A full charge allows up to 40 minutes of continuous use. **DO NOT** charge your *soniceraser pro* for more than 40 continuous hours. The Handle needs to receive a charge occasionally. If you rarely use the Handle, recharge it in the Charging Cradle at least every 30 days to maintain its charging capability. **NOTE:** Failure to regularly recharge the battery may cause the battery to stop recharging. We recommend that you keep the Handle in the Charging Cradle when not in use.
3. When using the *soniceraser pro* the Power Indicator Light will light up green until the battery reaches 20% remaining charge at which time the Power Indicator Light will light up red. When the remaining charge falls to 5% the Power Indicator Light will switch to a flashing red light.

## LIMITED WARRANTY

Please use this as your reference for any customer service or warranty needs you may have in the future.

**You must register your product within 90 days of purchase in order to obtain the warranty. Your warranty is not effective unless you register your product within the time provided. In order to register your product go on-line to [www.warrantymichaeltoddbeauty.com](http://www.warrantymichaeltoddbeauty.com) and complete the information required for registration.**

LIMITED 2 YEARS WARRANTY. This limited warranty gives you specific legal rights and you may have other rights, which vary from jurisdiction to jurisdiction. Any claim for breach of this limited warranty must be commenced within 2 years following the original purchase.

For a period of two(2) years from the date of purchase, Michael Todd Beauty ("MICHAEL TODD BEAUTY") warrants the *soniceraser pro* (the "Product") against Covered Defects. If there is a Covered Defect during the warranty period MICHAEL TODD BEAUTY agrees to exchange the Product with a new equivalent Product. This is the sole remedy available for breach of this limited warranty.

The limited warranty begins on the date of purchase and expires 2 years thereafter. In no event is the warranty period tolled or extended - including any period during which the Product (or any parts thereof) is in repair or if the Product is replaced by MICHAEL TODD BEAUTY during the warranty period.

This warranty is limited and applies only to a purchase of the Product in the United States from MICHAEL TODD BEAUTY or an authorized reseller, is valid for the original end user of the Product only and cannot be transferred under any circumstances. Except as otherwise provided by law, to obtain service under this warranty, return the Product with a Covered Defect to us with your purchase receipt. For shipping instructions, contact our Customer Service department [customerservice@michaeltoddbeauty.com](mailto:customerservice@michaeltoddbeauty.com) or call (772) 800-7654.

WHAT IS COVERED? "Covered Defects" means damage, defects or failure resulting from ordinary consumer use under normal home conditions according to the instruction manual included with the Product. Covered Defects does NOT include using the Product with a voltage converter; damage, defects or failure caused by misuse or improper care of the Product; impact or insertion of foreign objects into or on the Product; improper storage; improper or incorrectly performed maintenance or repair of the Product; application of harmful chemicals or abrasives to the Product; failure to follow manufacturer's instructions; normal wear and tear; alteration or overuse of the Product; as well as any professional or commercial use of the Product.

WARRANTY LIMITATIONS. THIS LIMITED WARRANTY IS EXCLUSIVE OF ALL OTHER WARRANTIES, WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED. THERE ARE NO IMPLIED WARRANTIES CREATED BY THE MANUFACTURE, SALE OR USE OF THE PRODUCT. IN NO EVENT SHALL MICHAEL TODD BEAUTY BE LIABLE FOR ANY SPECIAL, INCIDENTAL, CONSEQUENTIAL, PUNITIVE, OR INDIRECT DAMAGES BASED UPON BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, STRICT PRODUCT LIABILITY, TORT OR ANY OTHER LEGAL THEORY.

PRODUCT CHANGES. MICHAEL TODD BEAUTY reserves the right to make changes and improvements to its products without incurring any obligation to incorporate such improvements in products previously sold or to notify you of these changes or improvements.

DISPUTES. This limited warranty is governed by the laws of the State of Florida, without regard to principles of conflicts of laws, and applicable United States federal law. The United Nations Convention on Contracts for the International Sale of Goods shall not apply. YOU AGREE TO SUBMIT ANY CLAIM, DISPUTE, OR CONTROVERSY (REGARDLESS OF THE FORM OR TYPE OF ACTION) AGAINST US TO ARBITRATION ADMINISTERED BY THE NATIONAL ARBITRATION FORUM UNDER ITS CODE OF PROCEDURE THEN IN EFFECT (available at <http://www.arb-forum.com> or via telephone at 1-800-474-2371). Except to the extent otherwise required by law, you agree that the arbitration proceedings will take place exclusively where the principal office of MICHAEL TODD BEAUTY is located at the time of institution of the proceedings. You may be entitled to file suit under state laws without submitting your claim to arbitration first. Information on arbitration may be obtained from, and a claim may be filed with, the National Arbitration Forum at P.O. Box 5019, Minneapolis, MN 55405-0191 USA.

## TECHNICAL SPECIFICATIONS

Vibration Frequency: 36000 movements per minute

Temperature: 42° C

Unit Size: 7" with base 6.5" without base

Country of Origin: China

For more information, updates, reviews

frequently asked questions and more

visit us on the web at: [www.michaeltoddbeauty.com](http://www.michaeltoddbeauty.com)

This equipment has been tested and found to comply with  
Part 15 of the FCC Rules

Distributed by

Michael Todd Beauty LP

Port St Lucie, FL

[www.michaeltoddbeauty.com](http://www.michaeltoddbeauty.com) | 772.800.7654

*soniceraser pro* and Triplex Infusion Technology are  
trademarks of MICHAEL TODD BEAUTY LP.

© 2017 MICHAEL TODD BEAUTY

All rights reserved



**TRADE NAME: MICHAEL TODD BEAUTY**

**MODEL: SONICERASER PRO**

**Tested to Comply with FCC Standards**



